

For Research Use Only

BCDIN3D Recombinant antibody, PBS Only (Capture)

Catalog Number:83309-2-PBS



Basic Information

Catalog Number: 83309-2-PBS	GenBank Accession Number: BC053560	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 144233	CloneNo.: 240152E10
Source: Rabbit	UNIPROT ID: Q7Z5W3	
Isotype: IgG	Full Name: BCDIN3 domain containing	
Immunogen Catalog Number: AG15613	Calculated MW: 292 aa, 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00333-2, BCDIN3D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83309-2-PBS. Detection antibody: 83309-4-PBS. Standard: Ag15613. Range: 0.313-40 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83309-2-PBS against Human BCDIN3D was performed. The affinity constant is 0.864 nM.