

For Research Use Only

# SMARCA4 Recombinant antibody, PBS Only (Capture)



Catalog Number:83310-4-PBS

## Basic Information

Catalog Number:

83310-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16256

GenBank Accession Number:

BC150298

GeneID (NCBI):

6597

UNIPROT ID:

P51532

Full Name:

SWI/SNF related, matrix associated,  
actin dependent regulator of  
chromatin, subfamily a, member 4

Calculated MW:

1647 aa, 185 kDa

Purification Method:

Protein A purification

CloneNo.:

240172H1

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00334-3, SMARCA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83310-4-PBS. Detection antibody: 83310-3-PBS. Standard: Ag16256. Range: 0.625-40 ng/mL.