

For Research Use Only

BRAP Recombinant antibody, PBS Only (Detector)



Catalog Number: 83314-1-PBS

Basic Information

Catalog Number: 83314-1-PBS	GenBank Accession Number: BC136698	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8315	CloneNo.: 240099A2
Source: Rabbit	UNIPROT ID: Q7Z569	
Isotype: IgG	Full Name: BRCA1 associated protein	
Immunogen Catalog Number: AG26728	Calculated MW: 592 aa, 67 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00338-3, BRAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83314-3-PBS. Detection antibody: 83314-1-PBS. Standard: Ag26728. Range: 0.313-40 ng/mL