

For Research Use Only

BRAP Recombinant antibody, PBS Only



Catalog Number: **83314-6-PBS**

Basic Information

Catalog Number: 83314-6-PBS	GenBank Accession Number: BC136698	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8315	CloneNo.: 240099D8
Source: Rabbit	UNIPROT ID: Q7Z569	
Isotype: IgG	Full Name: BRCA1 associated protein	
Immunogen Catalog Number: AG26728	Calculated MW: 592 aa, 67 kDa	

Applications

Tested Applications:
IF, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed MCF-7 cells using BRAP antibody (83314-6-RR, Clone: 240099D8) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1). This data was developed using the same antibody clone with 83314-6-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83314-6-RR against Human BRAP were performed. The affinity constant is below 1 pM.