

For Research Use Only

NUPR1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83315-4-PBS

Basic Information

Catalog Number:

83315-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7039

GenBank Accession Number:

BC002434

GeneID (NCBI):

26471

UNIPROT ID:

O60356

Full Name:

nuclear protein 1

Calculated MW:

8 kDa

Purification Method:

Protein A purification

CloneNo.:

240149H8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00339-2, NUPR1 Recombinant Matched Antibody Pair - BSA and Azide Free. Capture antibody: 83315-2-PBS. Detection antibody: 83315-4-PBS. Standard: Ag7039. Range: 6.25-100 ng/mL.