

For Research Use Only

# ADRB2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83317-4-PBS

## Basic Information

Catalog Number:

83317-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3715

GenBank Accession Number:

BC012481

GeneID (NCBI):

154

UNIPROT ID:

P07550

Full Name:

adrenergic, beta-2-, receptor, surface

Calculated MW:

46 kDa

Purification Method:

Protein A purification

CloneNo.:

240186D9

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00341-3, ADRB2 Recombinant Matched Antibody Pair - BSA and Azide Free. Capture antibody: 83317-4-PBS. Detection antibody: 83317-3-PBS. Standard: Ag3715. Range: 0.156-20 ng/mL.