

For Research Use Only

APIP Recombinant antibody, PBS Only (Detector)



Catalog Number: 83320-2-PBS

Basic Information

Catalog Number: 83320-2-PBS	GenBank Accession Number: BC009077	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51074	CloneNo.: 240189H1
Source: Rabbit	UNIPROT ID: Q96GX9	
Isotype: IgG	Full Name: APAF1 interacting protein	
Immunogen Catalog Number: AG2729	Calculated MW: 242 aa, 27 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00346-3, AIP Recombinant Matched Antibody Pair - BSA and Azide Free. Capture antibody: 83320-3-PBS. Detection antibody: 83320-2-PBS. Standard: Ag2729. Range: 0.313-40 ng/mL.