

For Research Use Only

APIP Recombinant antibody, PBS Only (Capture)



Catalog Number: 83320-4-PBS

Basic Information

| | | |
|--|--|---|
| Catalog Number: 83320-4-PBS | GenBank Accession Number: BC009077 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 51074 | CloneNo.: 240189H7 |
| Source: Rabbit | UNIPROT ID: Q96GX9 | |
| Isotype: IgG | Full Name: APAF1 interacting protein | |
| Immunogen Catalog Number: AG2729 | Calculated MW: 242 aa, 27 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00346-2, AIP Recombinant Matched Antibody Pair - BSA and Azide Free. Capture antibody: 83320-4-PBS. Detection antibody: 83320-1-PBS. Standard: Ag2729. Range: 0.313-40 ng/mL.