

For Research Use Only

APIP Recombinant antibody

Catalog Number: 83320-5-RR



Basic Information

Catalog Number: 83320-5-RR	GenBank Accession Number: BC009077	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 51074	CloneNo.: 240189E2
Source: Rabbit	UNIPROT ID: Q96GX9	Recommended Dilutions: WB 1:5000-1:50000
Isotype: IgG	Full Name: APAF1 interacting protein	
Immunogen Catalog Number: AG2729	Calculated MW: 242 aa, 27 kDa	
	Observed MW: 30 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : PC-3 cells, K-562 cells.
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83320-5-RR (AP1P antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 83320-5-RR against Human AP1P were performed. The affinity constant is 5.86 nM.