

For Research Use Only

EGLN3 Recombinant antibody, PBS Only (Detector)



Catalog Number:83321-4-PBS

Basic Information

Catalog Number: 83321-4-PBS	GenBank Accession Number: BC010992	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 112399	CloneNo.: 240165A12
Source: Rabbit	UNIPROT ID: Q9H6Z9	
Isotype: IgG	Full Name: egl nine homolog 3 (C. elegans)	
Immunogen Catalog Number: AG27464	Calculated MW: 27 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00347-1, EGLN3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83321-1-PBS. Detection antibody: 83321-4-PBS. Standard: Ag27464. Range: 0.313-40 ng/mL.