

For Research Use Only

# NRG1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83323-1-PBS

## Basic Information

Catalog Number: 83323-1-PBS	GenBank Accession Number: BC007675	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3084	CloneNo.: 240184A4
Source: Rabbit	UNIPROT ID: Q02297	
Isotype: IgG	Full Name: neuregulin 1	
Immunogen Catalog Number: AG0803	Calculated MW: 70 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00232-3, NRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83323-4-PBS. Detection antibody: 83323-1-PBS. Standard: Ag0803. Range: 0.313-40 ng/mL