

For Research Use Only

NRG1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83323-3-PBS

Basic Information

Catalog Number: 83323-3-PBS	GenBank Accession Number: BC007675	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3084	CloneNo.: 240184C3
Source: Rabbit	UNIPROT ID: Q02297	
Isotype: IgG	Full Name: neuregulin 1	
Immunogen Catalog Number: AG0803	Calculated MW: 70 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00232-4, NRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83323-4-PBS. Detection antibody: 83323-3-PBS. Standard: Ag0803. Range: 0.313-40 ng/mL