

For Research Use Only

# NRG1 Recombinant antibody, PBS Only (Capture)



Catalog Number:83323-4-PBS

## Basic Information

|  |  |   |
|--|--|---|
| <b>Catalog Number:</b><br>83323-4-PBS      | <b>GenBank Accession Number:</b><br>BC007675 | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1 mg/ml                    | <b>GeneID (NCBI):</b><br>3084                | <b>CloneNo.:</b><br>240184A9                          |
| <b>Source:</b><br>Rabbit                   | <b>UNIPROT ID:</b><br>Q02297                 |   |
| <b>Isotype:</b><br>IgG                     | <b>Full Name:</b><br>neuregulin 1            |   |
| <b>Immunogen Catalog Number:</b><br>AG0803 | <b>Calculated MW:</b><br>70 kDa              |   |

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00232-3, NRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83323-4-PBS. Detection antibody: 83323-1-PBS. Standard: Ag0803. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00232-4, NRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83323-4-PBS. Detection antibody: 83323-3-PBS. Standard: Ag0803. Range: 0.313-40 ng/mL