

For Research Use Only

NRG1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83323-5-PBS

Basic Information

Catalog Number:

83323-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0803

GenBank Accession Number:

BC007675

GeneID (NCBI):

3084

UNIPROT ID:

Q02297

Full Name:

neuregulin 1

Calculated MW:

70 kDa

Purification Method:

Protein A purification

CloneNo.:

240184C8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00232-5, NRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83323-2-PBS. Detection antibody: 83323-5-PBS. Standard: Ag0803. Range: 0.313-40 ng/mL