

For Research Use Only

GANAB Recombinant antibody, PBS Only (Detector)

Catalog Number:83324-2-PBS



Basic Information

Catalog Number: 83324-2-PBS	GenBank Accession Number: BC017435	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23193	CloneNo.: 240279H5
Source: Rabbit	UNIPROT ID: Q14697	
Isotype: IgG	Full Name: glucosidase, alpha; neutral AB	
Immunogen Catalog Number: AG30417	Calculated MW: 944aa,107 kDa; 488aa,55 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00350-1, GANAB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83324-1-PBS. Detection antibody: 83324-2-PBS. Standard: Ag30417. Range: 1.25-40 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83324-2-PBS against Human GANAB was performed. The affinity constant is 0.653 nM.