

For Research Use Only

RAB13 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83336-1-PBS

Basic Information

Catalog Number: 83336-1-PBS	GenBank Accession Number: BC000799	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5872	CloneNo.: 240216C4
Source: Rabbit	UNIPROT ID: P51153	
Isotype: IgG	Full Name: RAB13, member RAS oncogene family	
Immunogen Catalog Number: AG34898	Calculated MW: 203 aa, 23 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00354-1, RAB13 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83336-2-PBS. Detection antibody: 83336-1-PBS. Standard: Ag34898. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00354-2, RAB13 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83336-3-PBS. Detection antibody: 83336-1-PBS. Standard: Ag34898. Range: 1.25-20 ng/mL