

For Research Use Only

DNASE1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83337-5-PBS

Basic Information

Catalog Number:

83337-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10117

GenBank Accession Number:

BC029437

GeneID (NCBI):

1773

UNIPROT ID:

P24855

Full Name:

deoxyribonuclease I

Calculated MW:

282 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

240222B11

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00356-1, DNASE1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83337-1-PBS. Detection antibody: 83337-5-PBS. Standard: Ag10117. Range: 0.313-40 ng/mL.