

For Research Use Only

ST3GAL5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83339-3-PBS

Basic Information

Catalog Number: 83339-3-PBS	GenBank Accession Number: BC065936	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8869	CloneNo.: 240148C10
Source: Rabbit	UNIPROT ID: Q9UNP4	
Isotype: IgG	Full Name: ST3 beta-galactoside alpha-2,3-sialyltransferase 5	
Immunogen Catalog Number: AG6207	Calculated MW: 48 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00381-1, ST3GALS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83339-2-PBS. Detection antibody: 83339-3-PBS. Standard: Ag6207. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00381-2, ST3GALS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83339-1-PBS. Detection antibody: 83339-3-PBS. Standard: Ag6207. Range: 0.625-40 ng/mL.