

For Research Use Only

# CRIP1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83340-2-PBS

## Basic Information

|  |   |   |
|--|---|---|
| <b>Catalog Number:</b><br>83340-2-PBS      | <b>GenBank Accession Number:</b><br>BC002738              | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1 mg/ml                    | <b>GeneID (NCBI):</b><br>1396                             | <b>CloneNo.:</b><br>240176H7                          |
| <b>Source:</b><br>Rabbit                   | <b>UNIPROT ID:</b><br>P50238                              |   |
| <b>Isotype:</b><br>IgG                     | <b>Full Name:</b><br>cysteine-rich protein 1 (intestinal) |   |
| <b>Immunogen Catalog Number:</b><br>AG7588 | <b>Calculated MW:</b><br>9 kDa                            |   |

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00382-1, CRP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83340-3-PBS. Detection antibody: 83340-2-PBS. Standard: Ag7588. Range: 0.156-20 ng/mL