

For Research Use Only

# MTF1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83344-2-PBS

## Basic Information

Catalog Number: 83344-2-PBS	GenBank Accession Number: BC014454	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4520	CloneNo.: 240336H1
Source: Rabbit	UNIPROT ID: Q14872	
Isotype: IgG	Full Name: metal-regulatory transcription factor 1	
Immunogen Catalog Number: AG21848	Calculated MW: 753 aa, 81 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00386-2, MTF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83344-2-PBS. Detection antibody: 83344-3-PBS. Standard: Ag21848. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP00386-3, MTF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83344-2-PBS. Detection antibody: 83344-1-PBS. Standard: Ag21848. Range: 0.625-20 ng/mL