

For Research Use Only

TSHR Recombinant antibody, PBS Only (Capture)



Catalog Number: 83348-2-PBS

Basic Information

Catalog Number: 83348-2-PBS	GenBank Accession Number: BC063613	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7253	CloneNo.: 240113G9
Source: Rabbit	UNIPROT ID: P16473	
Isotype: IgG	Full Name: thyroid stimulating hormone receptor	
Immunogen Catalog Number: AG5957	Calculated MW: 87 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00358-1, TSHR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83348-2-PBS. Detection antibody: 83348-3-PBS. Standard: Ag5957. Range: 0.625-80 ng/mL