

For Research Use Only

AGPAT2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83349-1-PBS

Basic Information

| | | |
|---|---|---|
| Catalog Number: 83349-1-PBS | GenBank Accession Number: BC019292 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 10555 | CloneNo.: 240153B6 |
| Source: Rabbit | UNIPROT ID: O15120 | |
| Isotype: IgG | Full Name: 1-acylglycerol-3-phosphate O-acyltransferase 2 (lysophosphatidic acid acyltransferase, beta) | |
| Immunogen Catalog Number: AG34761 | Calculated MW: 31 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00360-2, AGPAT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83349-3-PBS. Detection antibody: 83349-1-PBS. Standard: Ag34761. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00360-3, AGPAT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83349-4-PBS. Detection antibody: 83349-1-PBS. Standard: Ag34761. Range: 0.625-80 ng/mL.