

For Research Use Only

AGPAT2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83349-2-PBS

Basic Information

Catalog Number: 83349-2-PBS	GenBank Accession Number: BC019292	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10555	CloneNo.: 240153F11
Source: Rabbit	UNIPROT ID: O15120	
Isotype: IgG	Full Name: 1-acylglycerol-3-phosphate O-acyltransferase 2 (lysophosphatidic acid acyltransferase, beta)	
Immunogen Catalog Number: AG34761	Calculated MW: 31 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00360-1, AGPAT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83349-2-PBS. Detection antibody: 83349-4-PBS. Standard: Ag34761. Range: 0.625-80 ng/mL.