

For Research Use Only

CACHD1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83350-1-PBS

Basic Information

| | | |
|---|--|---|
| Catalog Number: 83350-1-PBS | GenBank Accession Number: BC039301 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 57685 | CloneNo.: 240197B2 |
| Source: Rabbit | UNIPROT ID: Q5VU97 | |
| Isotype: IgG | Full Name: cache domain containing 1 | |
| Immunogen Catalog Number: AG22138 | Calculated MW: 1274 aa, 142 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00362-1, CACHD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83350-2-PBS. Detection antibody: 83350-1-PBS. Standard: Ag22138. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00362-2, CACHD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83350-3-PBS. Detection antibody: 83350-1-PBS. Standard: Ag22138. Range: 0.625-80 ng/mL.