

For Research Use Only

CACHD1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83350-2-PBS

Basic Information

Catalog Number: 83350-2-PBS	GenBank Accession Number: BC039301	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57685	CloneNo.: 240197G8
Source: Rabbit	UNIPROT ID: Q5VU97	
Isotype: IgG	Full Name: cache domain containing 1	
Immunogen Catalog Number: AG22138	Calculated MW: 1274 aa, 142 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00362-1, CACHD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83350-2-PBS. Detection antibody: 83350-1-PBS. Standard: Ag22138. Range: 0.625-80 ng/mL.