

For Research Use Only

TCIRG1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83351-2-PBS

Basic Information

Catalog Number: 83351-2-PBS	GenBank Accession Number: BC018133	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10312	CloneNo.: 240166G12
Source: Rabbit	UNIPROT ID: Q13488	
Isotype: IgG	Full Name: T-cell, immune regulator 1, ATPase, H ⁺ transporting, lysosomal VO subunit A3	
Immunogen Catalog Number: AG3378	Calculated MW: 830 aa, 93 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00365-3, TCIRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83351-1-PBS. Detection antibody: 83351-2-PBS. Standard: Ag3378. Range: 0.313-40 ng/mL