

For Research Use Only

TCIRG1 Recombinant antibody, PBS Only



Catalog Number: 83351-6-PBS

Basic Information

Catalog Number: 83351-6-PBS	GenBank Accession Number: BC018133	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10312	CloneNo.: 240166B2
Source: Rabbit	UNIPROT ID: Q13488	
Isotype: IgG	Full Name: T-cell, immune regulator 1, ATPase, H+ transporting, lysosomal VO subunit A3	
Immunogen Catalog Number: AG3378	Calculated MW: 830 aa, 93 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83351-6-RR (TCIRG1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83351-6-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83351-6-RR against Human TCIRG1 were performed. The affinity constant is 2.14 nM.