

For Research Use Only

FGR Recombinant antibody, PBS Only (Detector)

Catalog Number: 83352-2-PBS



Basic Information

Catalog Number:

83352-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34944

GenBank Accession Number:

NM_005248

GeneID (NCBI):

2268

UNIPROT ID:

P09769

Full Name:

Gardner-Rasheed feline sarcoma viral (v-fgr) oncogene homolog

Calculated MW:

59 kDa

Purification Method:

Protein A purification

CloneNo.:

240275G6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00370-2, FGR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83352-3-PBS. Detection antibody: 83352-2-PBS. Standard: Ag34944. Range: 0.313-40 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83352-2-PBS against Human FGR was performed. The affinity constant is 0.291 nM.