

For Research Use Only

FGR Recombinant antibody, PBS Only (Capture)

Catalog Number:83352-4-PBS



Basic Information

Catalog Number: 83352-4-PBS	GenBank Accession Number: NM_005248	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2268	CloneNo.: 240275H1
Source: Rabbit	UNIPROT ID: P09769	
Isotype: IgG	Full Name: Gardner-Rasheed feline sarcoma viral (v-fgr) oncogene homolog	
Immunogen Catalog Number: AG34944	Calculated MW: 59 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00370-3, FGR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83352-4-PBS. Detection antibody: 83352-3-PBS. Standard: Ag34944. Range: 0.625-40 ng/mL



Bilayer interferometry (BLI) kinetic assay of 83352-4-PBS against Human FGR was performed. The affinity constant is below 1 pM.