

For Research Use Only

CNKS2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83354-1-PBS

Basic Information

Catalog Number: 83354-1-PBS	GenBank Accession Number: BC126121	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 22866	CloneNo.: 240350C8
Source: Rabbit	UNIPROT ID: Q8WXI2	
Isotype: IgG	Full Name: connector enhancer of kinase suppressor of Ras 2	
Immunogen Catalog Number: AG21174	Calculated MW: 1034 aa, 118 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00389-2, CNKSR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83354-1-PBS. Detection antibody: 83354-4-PBS. Standard: Ag21174. Range: 6.25-100ng/mL.