

For Research Use Only

SLC22A8 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83357-4-PBS

Basic Information

Catalog Number:

83357-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9527

GenBank Accession Number:

bc022387

GeneID (NCBI):

9376

UNIPROT ID:

Q8TCC7

Full Name:

solute carrier family 22 (organic anion transporter), member 8

Calculated MW:

542 aa, 60 kDa

Purification Method:

Protein A purification

CloneNo.:

240138G12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00363-1, SLC22A8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83357-4-PBS. Detection antibody: 83357-5-PBS. Standard: Ag9527. Range: 0.313-40ng/mL.