

For Research Use Only

SLC39A10 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83358-2-PBS

Basic Information

Catalog Number: 83358-2-PBS	GenBank Accession Number: BC101516	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57181	CloneNo.: 240261G3
Source: Rabbit	UNIPROT ID: Q9ULF5	
Isotype: IgG	Full Name: solute carrier family 39 (zinc transporter), member 10	
Immunogen Catalog Number: AG20099	Calculated MW: 831 aa, 94 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00367-1, SLC39A10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83358-3-PBS. Detection antibody: 83358-2-PBS. Standard: Ag20099. Range: 0.625-80ng/mL.



Cytometric bead array standard curve of MP00367-2, SLC39A10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83358-1-PBS. Detection antibody: 83358-2-PBS. Standard: Ag20099. Range: 0.625-80ng/mL.