

For Research Use Only

SLC25A46 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83360-3-PBS

Basic Information

Catalog Number: 83360-3-PBS	GenBank Accession Number: BC017169	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 91137	CloneNo.: 240114A9
Source: Rabbit	UNIPROT ID: Q96AG3	
Isotype: IgG	Full Name: solute carrier family 25, member 46	
Immunogen Catalog Number: AG27293	Calculated MW: 418 aa, 46 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00369-1, SLC25A46 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83360-2-PBS. Detection antibody: 83360-3-PBS. Standard: Ag27293. Range: 0.625-80ng/mL



Cytometric bead array standard curve of MP00369-2, SLC25A46 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83360-1-PBS. Detection antibody: 83360-3-PBS. Standard: Ag27293. Range: 2.5-80ng/mL