

For Research Use Only

KIAA1217 Recombinant antibody, PBS Only (Detector)



Catalog Number:83361-1-PBS

Basic Information

Catalog Number: 83361-1-PBS	GenBank Accession Number: BC018764	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 56243	CloneNo.: 240276A4
Source: Rabbit	UNIPROT ID: Q5T5P2	
Isotype: IgG	Full Name: KIAA1217	
Immunogen Catalog Number: AG19367	Calculated MW: 1943 aa, 214 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00371-1, KIAA1217 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83361-2-PBS. Detection antibody: 83361-1-PBS. Standard: Ag19367. Range: 0.78-100ng/mL.