

For Research Use Only

NCAM1 Recombinant antibody, PBS Only (Capture)



Catalog Number:83365-2-PBS

Basic Information

Catalog Number: 83365-2-PBS	GenBank Accession Number: BC047244	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4684	CloneNo.: 240078C12
Source: Rabbit	ENSEMBL Gene ID: ENSG00000149294	
Isotype: IgG	UNIPROT ID: P13591	
	Full Name: neural cell adhesion molecule 1	
	Calculated MW: 95 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00393-3, NCAM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83365-2-PBS. Detection antibody: 83365-1-PBS. Standard: Eg31816. Range: 2.5-40ng/mL.