

For Research Use Only

# Mouse VEGFR-2/Kdr Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: **83366-2-PBS**

## Basic Information

Catalog Number: 83366-2-PBS	GenBank Accession Number: NM_010612.2	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 16542	CloneNo.: 240214G8
Source: Rabbit	UNIPROT ID: P35918-1	
Isotype: IgG	Full Name: kinase insert domain protein receptor	
	Calculated MW: 150KD	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Mouse

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00394-2, MOUSE VEGFR2/KDR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83366-2-PBS. Detection antibody: 83366-3-PBS. Standard: Eg0703. Range: 2.5-80 ng/mL.



Cytometric bead array standard curve of MP00394-3, MOUSE VEGFR2/KDR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83366-1-PBS. Detection antibody: 83366-2-PBS. Standard: Eg0703. Range: 5-80 ng/mL.