

For Research Use Only

Mouse VEGFR-2/Kdr Recombinant antibody, PBS Only (Detector)



Catalog Number:83366-3-PBS

Basic Information

Catalog Number:

83366-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_010612.2

GeneID (NCBI):

16542

UNIPROT ID:

P35918-1

Full Name:

kinase insert domain protein receptor

Calculated MW:

150KD

Purification Method:

Protein A purification

CloneNo.:

240214F12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00394-1, MOUSE VEGFR2/KDR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83366-1-PBS. Detection antibody: 83366-3-PBS. Standard: Eg0703. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00394-2, MOUSE VEGFR2/KDR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83366-2-PBS. Detection antibody: 83366-3-PBS. Standard: Eg0703. Range: 2.5-80 ng/mL