

For Research Use Only

LIMS1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83368-1-PBS

Basic Information

Catalog Number: 83368-1-PBS	GenBank Accession Number: BC005341	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3987	CloneNo.: 240323C2
Source: Rabbit	UNIPROT ID: P48059	
Isotype: IgG	Full Name: LIM and senescent cell antigen-like domains 1	
Immunogen Catalog Number: AG14764	Calculated MW: 325 aa, 37 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00396-1, LIMS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83368-1-PBS. Detection antibody: 83368-2-PBS. Standard: Ag14764. Range: 1.25-80 ng/mL