

For Research Use Only

AHR Recombinant antibody, PBS Only (Capture)



Catalog Number: 83369-2-PBS

Basic Information

Catalog Number: 83369-2-PBS	GenBank Accession Number: BC070080	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 196	CloneNo.: 240325E9
Source: Rabbit	UNIPROT ID: P35869	
Isotype: IgG	Full Name: aryl hydrocarbon receptor	
Immunogen Catalog Number: AG12193	Calculated MW: 848 aa, 96 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00397-2, AHR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83369-2-PBS. Detection antibody: 83369-3-PBS. Standard: Ag12193. Range: 1.25-80 ng/mL