

For Research Use Only

# Rat MCP-1 Recombinant antibody, PBS Only (Capture)



Catalog Number:83370-3-PBS

## Basic Information

Catalog Number:

83370-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_031530

GeneID (NCBI):

24770

UNIPROT ID:

P14844

Full Name:

chemokine (C-C motif) ligand 2

Calculated MW:

16 kDa

Purification Method:

Protein A purification

CloneNo.:

240341E10

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Rat

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00390-1, RAT MCP-1/CCL-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83370-3-PBS. Detection antibody: 83370-1-PBS. Standard: Eg0397. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP00390-2, RAT MCP-1/CCL-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83370-3-PBS. Detection antibody: 83370-2-PBS. Standard: Eg0397. Range: 0.156-20 ng/mL