

For Research Use Only

BAP1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83371-2-PBS

Basic Information

Catalog Number: 83371-2-PBS	GenBank Accession Number: BC001596	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8314	CloneNo.: 240349F5
Source: Rabbit	UNIPROT ID: Q92560	
Isotype: IgG	Full Name: BRCA1 associated protein-1 (ubiquitin carboxy-terminal hydrolase)	
Immunogen Catalog Number: AG34415	Calculated MW: 81-91 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00357-1, BAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83371-1-PBS. Detection antibody: 83371-2-PBS. Standard: Ag34415. Range: 1.25-80 ng/mL