

For Research Use Only

KIAA1199 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83373-2-PBS

Basic Information

| | | |
|--------------------------------------|---------------------------------------|--|
| Catalog Number: 83373-2-PBS | GenBank Accession Number: BC020256 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 57214 | CloneNo.: 240372F3 |
| Source: Rabbit | UNIPROT ID: Q8WUJ3 | |
| Isotype: IgG | Full Name: KIAA1199 | |
| Immunogen Catalog Number: AG15527 | Calculated MW: 1361 aa, 153 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00373-2, KIAA1199 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83373-1-PBS. Detection antibody: 83373-2-PBS. Standard: Ag15527. Range: 1.25-80 ng/mL.