

For Research Use Only

Mouse TIGIT Recombinant antibody, PBS Only (Detector)



Catalog Number:83374-1-PBS

Basic Information

Catalog Number: 83374-1-PBS	GenBank Accession Number: NP_001139797.1	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 100043314	CloneNo.: 240242D6
Source: Rabbit	UNIPROT ID: A0A0B4J1G6	
Isotype: IgG	Full Name: T cell immunoreceptor with Ig and ITIM domains	
	Calculated MW: 26kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00374-2, MOUSE TIGIT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83374-2-PBS. Detection antibody: 83374-1-PBS. Standard: Eg0925. Range: 0.313-40 ng/mL.