

For Research Use Only

Mouse TIGIT Recombinant antibody, PBS Only (Capture)



Catalog Number: 83374-2-PBS

Basic Information

Catalog Number:

83374-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NP_001139797.1

GeneID (NCBI):

100043314

UNIPROT ID:

A0A0B4J1G6

Full Name:

T cell immunoreceptor with Ig and

ITIM domains

Calculated MW:

26kDa

Purification Method:

Protein A purification

CloneNo.:

240242G1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00374-1, MOUSE TIGIT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83374-2-PBS. Detection antibody: 83374-3-PBS. Standard: Eg0925. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00374-2, MOUSE TIGIT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83374-2-PBS. Detection antibody: 83374-1-PBS. Standard: Eg0925. Range: 0.313-40 ng/mL.