

For Research Use Only

NEIL1 Recombinant antibody, PBS Only (Capture)



Catalog Number:83376-4-PBS

Basic Information

Catalog Number: 83376-4-PBS	GenBank Accession Number: BC010876	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79661	CloneNo.: 240322E8
Source: Rabbit	UNIPROT ID: Q96FI4	
Isotype: IgG	Full Name: nei endonuclease VIII-like 1 (E. coli)	
Immunogen Catalog Number: AG2793	Calculated MW: 390 aa, 44 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00380-2, NEIL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83376-4-PBS. Detection antibody: 83376-3-PBS. Standard: Ag2793. Range: 0.78-100 ng/mL