

For Research Use Only

RPL29 Recombinant antibody, PBS Only



Catalog Number: 83377-2-PBS

Basic Information

Catalog Number: 83377-2-PBS	GenBank Accession Number: BC008926	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6159	CloneNo.: 240062A3
Source: Rabbit	UNIPROT ID: P47914	
Isotype: IgG	Full Name: ribosomal protein L29	
Immunogen Catalog Number: AG34100	Calculated MW: 159 aa, 18 kDa	

Applications

Tested Applications:
WB, IF, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83377-2-RR (RPL29 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83377-2-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using RPL29 antibody (83377-2-RR, Clone: 240062A3) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 83377-2-PBS in a different storage buffer formulation.