

For Research Use Only

Cytokeratin 14 Recombinant antibody, PBS Only



Catalog Number:83379-1-PBS

Basic Information

Catalog Number: 83379-1-PBS	GenBank Accession Number: BC002690	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3861	CloneNo.: 1H22
Source: Rabbit	UNIPROT ID: P02533	
Isotype: IgG	Full Name: keratin 14	
Immunogen Catalog Number: AG0188	Calculated MW: 472 aa, 52 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat skin tissue were subjected to SDS PAGE followed by western blot with 83379-1-RR (KRT14 antibody) at dilution of 1:1400 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83379-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83379-1-RR against Human Cytokeratin 14 were performed. The affinity constant is below 1 pM.