

For Research Use Only

CD28 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83383-3-PBS

Basic Information

Catalog Number:

83383-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC093698

GeneID (NCBI):

940

ENSEMBL Gene ID:

ENSG00000178562

UNIPROT ID:

P10747

Full Name:

CD28 molecule

Calculated MW:

220 aa, 25 kDa

Purification Method:

Protein A purification

CloneNo.:

240354E4

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00136-4, CD28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83383-3-PBS. Detection antibody: 83383-2-PBS. Standard: Eg0889. Range: 0.313-40 ng/mL.