

For Research Use Only

CD28 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83383-4-PBS

Basic Information

Catalog Number: 83383-4-PBS	GenBank Accession Number: BC093698	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 940	CloneNo.: 240354F2
Source: Rabbit	ENSEMBL Gene ID: ENSG00000178562	
Isotype: IgG	UNIPROT ID: P10747	
	Full Name: CD28 molecule	
	Calculated MW: 220 aa, 25 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00136-3, CD28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83383-1-PBS. Detection antibody: 83383-4-PBS. Standard: Eg0889. Range: 0.313-40 ng/mL.