

For Research Use Only

# NLRC5 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83384-3-PBS

## Basic Information

<b>Catalog Number:</b> 83384-3-PBS	<b>GenBank Accession Number:</b> BC063566	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 84166	<b>CloneNo.:</b> 240338B12
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q86WI3	
<b>Isotype:</b> IgG	<b>Full Name:</b> NLR family, CARD domain containing 5	
<b>Immunogen Catalog Number:</b> AG34508	<b>Calculated MW:</b> 204 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



Cytometric bead array standard curve of MP00402-1, NLRC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83384-3-PBS. Detection antibody: 83384-4-PBS. Standard: Ag34508. Range: 1.56-100 ng/mL



Cytometric bead array standard curve of MP00402-2, NLRC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83384-1-PBS. Detection antibody: 83384-3-PBS. Standard: Ag34508. Range: 1.56-100 ng/mL



Cytometric bead array standard curve of MP00402-3, NLRC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83384-2-PBS. Detection antibody: 83384-3-PBS. Standard: Ag34508. Range: 1.56-100 ng/mL